# FlowCell

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**HAMILTON** 

# Table of Contents

<b>1</b>	<b>Introduction</b>	<b>4</b>
<b>2</b>	<b>Safety instructions</b>	<b>4</b>
2.1	General Safety instructions	4
2.2	Other risks	5
2.3	Intended use	5
2.4	Basic safety hints	5
<b>3</b>	<b>Transportation and storage</b>	<b>6</b>
<b>4</b>	<b>Description of the device</b>	<b>6</b>
<b>5</b>	<b>Assembly</b>	<b>8</b>
5.1	General	8
5.2	Mounting the FlowCell	8
5.3	Assembling the FlowCell	8
<b>6</b>	<b>Operation</b>	<b>9</b>
<b>7</b>	<b>Cleaning and maintenance</b>	<b>9</b>
<b>8</b>	<b>Technical data</b>	<b>11</b>
<b>9</b>	<b>Ordering information</b>	<b>12</b>

## Important note

Copyright © 2016 Hamilton Bonaduz AG, Bonaduz Switzerland. All rights reserved. All other trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

The reproduction of any part of this document in any form is forbidden without the express written agreement of Hamilton Bonaduz AG.

Contents of this manual can be modified without previous announcement. Technical modifications reserved. Greatest possible care was used on the correctness of the information in this manual. If errors should be discovered nevertheless, Hamilton Bonaduz AG is pleased to be informed about it. Regardless of this, Hamilton Bonaduz AG cannot assume liability for any errors in this manual or for their consequences.

# 1 Introduction

These operating instructions are intended for the FlowCell in the following versions:

Name	Ref
FlowCell TC25	242585-XYZ
FlowCell TC50	242590-XYZ

Hamilton flow-through cells are quality products, produced in accordance with the latest scientific and technical findings. Following the instructions given here assures maximum safety and product longevity.

These instructions must be read, understood and followed by all staff using the device. Hamilton assumes no responsibility for damage and for operational disruptions arising from failure to observe these instructions.

## 2 Safety Instructions

### 2.1 General Safety instructions

This flow-through cell and its parts are built using state-of-the-art technology and safe to operate. However, there can be risk to life and limb if users do not operate them correctly and appropriately.

The flow-through cell must only be used:

- For its intended purpose
- In optimum safety and running condition

Everyone involved in the operation and maintenance of this flow-through cell must carefully read and follow these operating instructions. Your personal safety may be at stake.

The process pressure of 16 bar and the process temperature of 140 °C must not be exceeded in any circumstances.

### 2.2 Other risks

Even after all necessary safety measures have been complied with, a risk of leaking or mechanical damage remains to the FlowCell. In any place where there are seals, there is a potential for gases or liquids to leak out undetectedly.

### 2.3 Intended use

The FlowCell is designed for specific operating conditions (see chapter 9) with respect to temperature, pressure and mounting.

The intended use also includes conformity with the conditions set down by Hamilton for assembly, initial operation, operation and maintenance. Any other or further use is classified as unintended, and Hamilton accepts no responsibility for any ensuing damage.

### 2.4 Basic safety hints

Installation, initial operation, cleaning, maintenance, troubleshooting

- Prescribed maintenance and inspection tasks must be completed at appropriate intervals, given in the chapter “Cleaning and Maintenance”.
- Operators should be informed before the start about maintenance and servicing.
- All attachments and media used with the FlowCell must be protected against unsupervised use.
- Sensors in the FlowCell must be dismantled only at ambient pressure.
- Defective parts must be replaced immediately. Use only original spare parts. There is no guarantee that parts from other manufacturers meet the operational and safety requirements of the FlowCell.
- On completing maintenance, check that the inside of the FlowCell is clean, and all the screw connections are screwed in tight.

### 3 Transportation and storage

Always do the following:

- Check the FlowCell for completeness immediately on receipt.
- Notify in writing without delay any damage in transit or any missing part.
- Before returning the FlowCell, contact Hamilton for a Return of Good Authorisation (RGA) number. Please make sure that the FlowCell is cleaned (decontaminated) before returning.
- When returning a FlowCell to a Hamilton partner or directly to Hamilton, make sure it is in its original packaging.
- When storing, protect the FlowCell from wetness.

### 4 Description of the FlowCell

The FlowCell is intended for the professional use in process analytics, as an adapter to mount sensors with a shaft diameter of 12 mm. It can be mounted into a pipe by means of Triclamp or Swagelok connectors.

Compatibility Table for Sensors:

	Tube A	Tube B
<b>Position 1</b>	pH Sensor	Conducell 4USF
	Conducell UPW *	OxyFerm
		VisiFerm
		VisiTrace
		VisiPro
<b>Position 2</b>		Conducell 4USF
		OxyFerm
		VisiFerm
		VisiTrace
		VisiPro

*\* not for Flowcell 242585 with Swagelok connection*


Figure 1: Components of the FlowCell 242585 with two sensor tubes

Figure 2: Components of the FlowCell 242585 with one sensor tube


Figure 3: Components of the FlowCell 242590 with two sensor tubes

Figure 4: Components of the FlowCell 242590 with one sensor tube

## 5 Assembly

### 5.1 General

Hamilton accepts no responsibility for damage resulting from incorrect assembly work by third parties.

If a FlowCell is modified, Hamilton does not accept liability and does not provide any further guarantee for the FlowCell.

### 5.2 Mounting the FlowCell

The FlowCell should always be mounted in horizontal position (see picture). It must not be mounted in vertical position, because this can prevent proper function of the FlowCell.

The body of the FlowCell can be easily mounted to the wall or a board with the help of 2 screws included with the FlowCell.

The FlowCell must be connected by means of Triclamp or Swagelok connectors. The flow should always lead from connection 1 to connection 2 (see picture). This way, a minimum volume of air and selfdraining of the FlowCell is ensured.

### 5.3 Assembling the Sensor

Ensure that there is no damage to the sensor or the FlowCell. Check whether all o-rings (see “Cleaning and Maintenance”) are in place in the appropriate grooves and are free of damage.

To avoid any mechanical damage to o-rings on assembly, they should be lightly greased.

When mounting two sensors, make certain they are mounted into the correct sensor tube (see compatibility table in chapter “Description of the FlowCell”).


Figure 5: Mounting of FlowCell

## 6 Operation

In EU countries, initial operation of the FlowCell ist permitted only when it has been verified that the FlowCell conforms to the requirement of the EU Directive.

Before initial operation of the FlowCell

- Ensure that the seal is tight
- Ensure that all parts are in working order.

In addition to initial operation, perform the above checks:

- After any breakdown
- When the device has been out of use for a long period
- After repairs or maintenance

**⚠ ATTENTION! These specifications must not be exceeded:  
Process temperature: -10°C to 140°C. Maximum process pressure: 16 bar.**

## 7 Cleaning and Maintenance

When carrying out all maintenance, cleaning and servicing, the general safety precautions given in this set of operating instructions must be observed. Operation safety and life expectancy of the flow-through cell are both dependent on appropriate maintenance and servicing, among other factors.

The FlowCell can be cleaned with standard cleaning agents. Never use scouring agents which contain hard particles. Please check the o-rings regularly. When mounting a new o-ring take care not to damage the o-ring nut or the o-ring itself. Otherwise the tightness can no longer be guaranteed.

Figure 6: Cut through FlowCell 242590.

All o-rings should be replaced regularly:

Figure 7: FlowCell with marked sensor positions and FlowCell disassembled without PEEK insert.

Figure 8: Swagelok connector enlarged.

NOTE: The pH sensor must always be mounted in position 1.

NOTE: An attempt to disassemble the Swagelok connector can lead to leakage of the FlowCell. Thus, the Swagelok connector must not be disassembled.

## 8 Technical Data

Process connection	Triclamp or Swagelok
Wetted parts	PEEK LSG / FG nature, Stainless Steel 1.4435
Non wetted parts	Stainless Steel 1.4305
Standard seals	EPDM (FDA approved)
Temperature range	-10 – 140 °C
Maximum pressure	0 – 16 bar
Internal volume Ref 242585	approx. 8 mL (only within the Peek cell)
Internal volume Ref 242590	approx. 25 mL (only within the Peek cell)
ATEX approval	Conform to DIN EN 13463-1

# 9 Ordering Information

## FlowCell

242585

Code	Measuring position
1	only Tube A (short)
2	only Tube B (long)
3	Tube A (short) and Tube B (long)
4	2 x Tube B (long)
0	special
Code	Pipe Connection
1	1/4" TC25
2	3/8" TC25
3	1/2" TC25
4	Swagelok 6 mm
5	Swagelok 10 mm *)
6	Swagelok 1/4"
7	Swagelok 3/8" *)
8	Swagelok 1/2" *)
0	special
Code	O-ring material
1	EPDM
2	FFPM (two measuring positions)
3	FFPM (one measuring position)
0	special

242585 -	← Ordercode
----------	-------------

\*) not self draining

### Ref Accessories

237387	O-ring kit FlowCell TC25 (EPDM)
242540	Sensor Dummy 96 (short)
242563	Sensor Dummy 117 (long)

Figure 9: FlowCell 242585 and connectors.

**FlowCell**

**242590**

Code	Measuring position
1	only Tube A (short)
2	only Tube B (long)
3	Tube A (short) and Tube B (long)
4	2 x Tube B (long)
0	special

Code	Pipe Connection
1	¾" TC50
2	1" TC50
3	1.5" TC50 *)
0	special

Code	o-ring material
1	EPDM
2	FFPM (two measuring positions)
3	FFPM (one measuring position)
0	special

<b>242590 -</b>	<b>← Ordercode</b>
-----------------	--------------------

\*) not self draining

**Ref Accessories**

<b>237390</b>	O-ring kit FlowCell TC50 (EPDM)
<b>242540</b>	Sensor Dummy 96 (short)
<b>242563</b>	Sensor Dummy 117 (long)

Figure 10: FlowCell 242590 and connectors.

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Lined area for notes on the left side of the page.

Lined area for notes on the right side of the page.

© 2016 Hamilton Bonaduz AG. All rights reserved.

REF 624521/01 — 03/2016

**HAMILTON**®

Web: [www.hamiltoncompany.com](http://www.hamiltoncompany.com)

USA: 800-648-5950

Europe: +41-58-610-10-10

**Hamilton Americas & Pacific Rim**

4970 Energy Way  
Reno, Nevada 89502 USA  
Tel: +1-775-858-3000  
Fax: +1-775-856-7259  
[sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

**Hamilton Europe, Asia & Africa**

Via Crusch 8  
CH-7402 Bonaduz, Switzerland  
Tel: +41-58-610-10-10  
Fax: +41-58-610-00-10  
[contact.pa.ch@hamilton.ch](mailto:contact.pa.ch@hamilton.ch)

To find a representative in your area, please visit [www.hamiltoncompany.com](http://www.hamiltoncompany.com).

This guide may be available in other languages. Visit [www.hamiltoncompany.com](http://www.hamiltoncompany.com) for more information.

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Copyright © 2016 Hamilton Bonaduz AG, Bonaduz Schweiz. Alle Rechte vorbehalten.