# Intelligent pH/ORP Transmitters

PC-3310/PC-3310-RS

21 3681 7199 - 21 99607 2513 contato@intertechrio.com.br www.intertechrio.com.br



L:Back △: ▲ L: ▲ EMT: Enter

Enter

### **Applications:**

- **Industrial Process**
- Chemical Process
- ♦ Water & Wastewater Treatment ▶ Environmental protection
- Scrubber
- Biochemistry

#### **Features:**

- 6 144 x 144 mm Intelligent pH/ORP Transmitter. €
- Waterproof and dustproof design
- 6 Solution-ground design eliminate interference of ground potential
- Large LCM with auto-sense backlight and contrast setting function
- 6 User-friendly operation with text and graphic illustration, clear signs and easy to understand
- 6 Multi-Cal calibration, Tech and NIST buffer automatic recognition
- Max. 3-point calibration, the calibration line calculated by linear regression, the last calibration data displayable
- With password and multi-key design-security protection enhancement
- Automatic temperature compensation, NTC30K or PT 1000 selectable
- Two sets of isolated current output DC 0/4~20mA, max. load 500Ω
- 6 Hi/Lo, Hi/Hi, Lo/Lo selectable two limited contact outputs with programmable set-points and hysteresis
- 6 Gradually continuous sample measurements average function to enhance stability
- ♦ Automatic sensor clean function with programmable time interval
- Three separate LED indicated alarms for recognition at a glance even from far place
- ♦ IP 65 (NEMA 4X) approved

#### **PC-3310-RS Only:**

- Measuring value can be selected for figure or chart display
- 6 Trace mode can display historical measured value curve from 3-minute to 4-week program
- Equipped with logbook function
- Max. 5-point calibration, with customer definition buffer setting and calibration curve, which can be displayed
- ♦ Equipped with one analog output and RS-485 interface
- Product adjustment function for field correction





## Intelligent pH/ORP Transmitters PC-3310/PC-3310-RS

#### **Specifications:**

Model		PC-3310	PC-3310-RS
Measuring mode		pH/ORP/Temp.	
Ranges	pН	-2.00~16.00 pH	
	ORP	-1999~1999 mV	
	Temp.	-30.0~130.0°C	
Resolution	pН	0.01 pH	
	ORP	1 mV	
	Temp.	0.1°C	
Accuracy	pН	± 0.01 pH ± 1 digit	
	ORP	± 0.1% ±1 digit	
	Temp.	± 0.2°C ±1 digit with temperature error correction function	
Temperature compensation		Automatic with NTC 30KΩ or PT-1000	
		Manual adjustment	
Calibration mode		Tech, NIST & Any, up to three points	Tech, NIST, Any & Define buffers, up to five points
Ambient Temperature		0~50°C	
Storage Temperature		-20~70°C	
Input Impedance		$>10^{12}\Omega$	
Display		Large LCM with auto-sense backlight and contrast function	
		Text mode	Text, Chart or Trace mode
Language		English, Traditional Chinese & Simplified Chinese	
Real-time chart mode		_	3 minutes real-time dynamic graph
Analog output 1		Isolated DC 0/4~20mA corresponding to pH/ORP, max. load 500 $\Omega$	
Analog output 2		Isolated DC 0/4~20mA corresponding to Temp., max. load $500\Omega$	_
Logbook		_	50 data records
Serial Interface		_	RS-485(MODBUS RTU or ASCII)
Setting	Contact	240VAC, 0.5A max. (recommend)	
	Activate	Hi/Lo, Hi/Hi, Lo/Lo selectable two limited programmable, ON/OFF	
Wash	Contact	240VAC, 0.5A max. (recommend)	
VV4311	Time	ON: 0~99min.59sec. OFF: 0~999 hours 59min.	
Voltage output		DC ± 12V, 1W max.	
Power Supply		100~240VAC±10%, 6W max., 50/60 Hz	
Installation		Wall or pipe or panel mounting	
Dimensions		144 mm x 144 mm x 115 mm (H x W x D)	
Cut off Dimensions		138 mm x 138 mm (H x W)	
Weight		0.8 kg	
Protection		IP 65 (NEMA 4X)	

<sup>\*</sup>The specifications and appearance of the instrument are subject to change without notice.

#### **Optional Accessories:**

1.pH electrode 9.Junction box 2.ORP electrode

10.Flow-through chamber

3.PT-1000 ATC probe 11.Immersion electrode holder 1M, 2M and 3M  $4.NTC~30K\Omega~ATC~probe$ 12.Industrial electrode protection housing/fitting

5. Field indication transmitter 13.Extension cable 6.pH buffer solutions 14.Sensor cleaning kit

7.ORP buffer solution 15.Sun shield

8. Cleaning solution 16.Stainless holder support base



Front clamshell design for easy wiring



